



DNASE1 Antibody, FITC conjugated

Product Code	CSB-PA007049LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24855
Immunogen	Recombinant Human Deoxyribonuclease-1 protein (78-231AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyribonuclease-1 (EC 3.1.21.1) (Deoxyribonuclease I) (DNase I) (Dornase alfa), DNASE1, DNL1 DRNI
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DNASE1