



RPS5 Antibody, FITC conjugated

Product Code	CSB-PA00729C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46782
Immunogen	Recombinant Human 40S ribosomal protein S5 protein (1-204AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S5 (Small ribosomal subunit protein uS7) [Cleaved into: 40S ribosomal protein S5, N-terminally processed], RPS5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPS5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.