



EGR2 Antibody, FITC conjugated

Product Code	CSB-PA007485LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11161
Immunogen	Recombinant Human E3 SUMO-protein ligase EGR2 protein (42-245AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 SUMO-protein ligase EGR2 (EC 2.3.2.-) (AT591) (E3 SUMO-protein transferase ERG2) (Early growth response protein 2) (EGR-2) (Zinc finger protein Krox-20), EGR2, KROX20
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EGR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.