



ESM1 Antibody, Biotin conjugated

Product Code	CSB-PA007825HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ30
Immunogen	Recombinant Human Endothelial cell-specific molecule 1 protein (20-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endothelial cell-specific molecule 1 (ESM-1), ESM1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ESM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.