



F12 Antibody

| | |
|---------------------|--|
| Product Code | CSB-PA007918LA01PI |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O97507 |
| Immunogen | Recombinant Pig Coagulation factor XII protein (20-371AA) |
| Raised In | Rabbit |
| Species Reactivity | Pig |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Coagulation factor XII (EC 3.4.21.38) (Hageman factor) (HAF) [Cleaved into: Coagulation factor XIIa heavy chain; Coagulation factor XIIa light chain], F12 |
| Immunogen Species | Sus scrofa (Pig) |
| Research Area | Others |
| Target Names | F12 |