



F12 Antibody, HRP conjugated

Product Code	CSB-PA007918LB01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O97507
Immunogen	Recombinant Pig Coagulation factor XII protein (20-371AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coagulation factor XII (EC 3.4.21.38) (Hageman factor) (HAF) [Cleaved into: Coagulation factor XIIa heavy chain; Coagulation factor XIIa light chain], F12
Immunogen Species	Sus scrofa (Pig)
Research Area	Others
Target Names	F12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.