



F2RL1 Antibody, Biotin conjugated

Product Code	CSB-PA007925LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55085
Immunogen	Recombinant Human Proteinase-activated receptor 2 protein (348-397AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteinase-activated receptor 2 (PAR-2) (Coagulation factor II receptor-like 1) (G-protein coupled receptor 11) (Thrombin receptor-like 1) [Cleaved into: Proteinase-activated receptor 2, alternate cleaved 1; Proteinase-activated receptor 2, alternate cleaved 2], F2RL1, GPR11 PAR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	F2RL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.