





FCGR3A Antibody

Product Code	CSB-PA008543LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08637
Immunogen	Recombinant Human Low affinity immunoglobulin gamma Fc region receptor III-A protein (49-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin gamma Fc region receptor III-A (CD16a antigen) (Fc-gamma RIII-alpha) (Fc-gamma RIII) (Fc-gamma RIIIa) (FcRIII) (FcRIIIa) (FcR-10) (IgG Fc receptor III-2) (CD antigen CD16a), FCGR3A, CD16A FCG3 FCGR3 IGFR3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCGR3A


Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA008543LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA008543LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA008543LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)