



FCGR3A Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA008543LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P08637 |
| Immunogen | Recombinant Human Low affinity immunoglobulin gamma Fc region receptor III-A protein (49-184AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Low affinity immunoglobulin gamma Fc region receptor III-A (CD16a antigen) (Fc-gamma RIII-alpha) (Fc-gamma RIII) (Fc-gamma RIIIa) (FcRIII) (FcRIIIa) (FcR-10) (IgG Fc receptor III-2) (CD antigen CD16a), FCGR3A, CD16A FCG3 FCGR3 IGFR3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | FCGR3A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |