



FCGR3A Antibody, Biotin conjugated

Product Code	CSB-PA008543LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08637
Immunogen	Recombinant Human Low affinity immunoglobulin gamma Fc region receptor III-A protein (49-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin gamma Fc region receptor III-A (CD16a antigen) (Fc-gamma RIII-alpha) (Fc-gamma RIII) (Fc-gamma RIIIa) (FcRIII) (FcRIIIa) (FcR-10) (IgG Fc receptor III-2) (CD antigen CD16a), FCGR3A, CD16A FCG3 FCGR3 IGFR3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCGR3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.