



FDPS Antibody, Biotin conjugated

Product Code	CSB-PA008563LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14324
Immunogen	Recombinant Human Farnesyl pyrophosphate synthase protein (1-419AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Farnesyl pyrophosphate synthase (FPP synthase) (FPS) (EC 2.5.1.10) ((2E,6E)-farnesyl diphosphate synthase) (Dimethylallyltranstransferase) (EC 2.5.1.1) (Farnesyl diphosphate synthase) (Geranyltranstransferase), FDPS, FPS KIAA1293
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FDPS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.