



fgf3 Antibody, HRP conjugated

Product Code	CSB-PA008630LB01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48802
Immunogen	Recombinant Zebrafish Fibroblast growth factor 3 protein (19-256AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 3 (FGF-3) (Heparin-binding growth factor 3) (HBGF-3), fgf3, fgf-3
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	fgf3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.