



FGG Antibody, FITC conjugated

Product Code	CSB-PA008651LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02679
Immunogen	Recombinant Human Fibrinogen gamma chain protein (278-437AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibrinogen gamma chain, FGG
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.