



FGR Antibody, FITC conjugated

Product Code	CSB-PA008655LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09769
Immunogen	Recombinant Human Tyrosine-protein kinase Fgr protein (11-72AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase Fgr (EC 2.7.10.2) (Gardner-Rasheed feline sarcoma viral (v-fgr) oncogene homolog) (Proto-oncogene c-Fgr) (p55-Fgr) (p58-Fgr) (p58c-Fgr), FGR, SRC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.