



FMOD Antibody, Biotin conjugated

Product Code	CSB-PA008755LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06828
Immunogen	Recombinant Human Fibromodulin protein (18-89AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibromodulin (FM) (Collagen-binding 59 kDa protein) (Keratan sulfate proteoglycan fibromodulin) (KSPG fibromodulin), FMOD, FM SLRR2E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FMOD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.