




FUCA1 Antibody

Product Code	CSB-PA009062LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04066
Immunogen	Recombinant Human Tissue alpha-L-fucosidase protein (32-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tissue alpha-L-fucosidase (EC 3.2.1.51) (Alpha-L-fucosidase I) (Alpha-L-fucoside fucosylhydrolase 1) (Alpha-L-fucosidase 1), FUCA1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FUCA1

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA009062LA01HU at dilution of 1:100