



FUCA1 Antibody, Biotin conjugated

Product Code	CSB-PA009062LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04066
Immunogen	Recombinant Human Tissue alpha-L-fucosidase protein (32-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tissue alpha-L-fucosidase (EC 3.2.1.51) (Alpha-L-fucosidase I) (Alpha-L-fucoside fucohydrolase 1) (Alpha-L-fucosidase 1), FUCA1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FUCA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.