




GALK1 Antibody

Product Code	CSB-PA009198LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51570
Immunogen	Recombinant Human Galactokinase protein (23-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galactokinase (EC 2.7.1.6) (Galactose kinase), GALK1, GALK
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GALK1

Image



Western Blot

Positive WB detected in: HEK293 whole cell lysate, HeLa whole cell lysate


All lanes: GALK1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 43, 46 kDa

Observed band size: 43 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA009198LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.