



GAS7 Antibody, HRP conjugated

Product Code	CSB-PA009271LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60861
Immunogen	Recombinant Human Growth arrest-specific protein 7 protein (1-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth arrest-specific protein 7 (GAS-7), GAS7, KIAA0394
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GAS7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.