




Gfap Antibody

Product Code	CSB-PA009369LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47819
Immunogen	Recombinant Rat Glial fibrillary acidic protein (1-430AA)
Raised In	Rabbit
Species Reactivity	Rat, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glial fibrillary acidic protein (GFAP), Gfap
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Gfap

Image



Western Blot

Positive WB detected in: Mouse brain tissue, Rat brain tissue

All lanes: Gfap antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50, 49 kDa

Observed band size: 50 kDa