



Gfap Antibody, HRP conjugated

Product Code	CSB-PA009369LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47819
Immunogen	Recombinant Rat Glial fibrillary acidic protein (1-430AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glial fibrillary acidic protein (GFAP), Gfap
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Gfap
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.