



Gfap Antibody, FITC conjugated

Product Code	CSB-PA009369LC01RA
Abbreviation	Glial fibrillary acidic protein
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47819
Immunogen	Recombinant Rat Glial fibrillary acidic protein (1-430AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Relevance	GFAP, a class-III intermediate filament, is a cell-specific marker that, during the development of the central nervous system, distinguishes astrocytes from other glial cells.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glial fibrillary acidic protein (GFAP), Gfap
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Gfap