



GGH Antibody, FITC conjugated

Product Code	CSB-PA009389LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92820
Immunogen	Recombinant Human Gamma-glutamyl hydrolase protein (219-293AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-glutamyl hydrolase (EC 3.4.19.9) (Conjugase) (GH) (Gamma-Glu-X carboxypeptidase), GGH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GGH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.