



GIPC1 Antibody, HRP conjugated

Product Code	CSB-PA009435LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14908
Immunogen	Recombinant Human PDZ domain-containing protein GIPC1 protein (7-74AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDZ domain-containing protein GIPC1 (GAIP C-terminus-interacting protein) (RGS-GAIP-interacting protein) (RGS19-interacting protein 1) (Synectin) (Tax interaction protein 2) (TIP-2), GIPC1, C19orf3 GIPC RGS19IP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GIPC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.