




# GIT1 Antibody

|                            |                                                                                                                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA009439LA01HU                                                                                                                                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                               |
| <b>Uniprot No.</b>         | Q9Y2X7                                                                                                                                                                                                      |
| <b>Immunogen</b>           | Recombinant Human ARF GTPase-activating protein GIT1 protein (470-642AA)                                                                                                                                    |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                      |
| <b>Species Reactivity</b>  | Human, Mouse                                                                                                                                                                                                |
| <b>Tested Applications</b> | ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200                                                                                                                    |
| <b>Form</b>                | Liquid                                                                                                                                                                                                      |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                                                                                              |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                                                            |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                                                    |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                  |
| <b>Alias</b>               | ARF GTPase-activating protein GIT1 (ARF GAP GIT1) (Cool-associated and tyrosine-phosphorylated protein 1) (CAT-1) (CAT1) (G protein-coupled receptor kinase-interactor 1) (GRK-interacting protein 1), GIT1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                        |
| <b>Research Area</b>       | Signal Transduction                                                                                                                                                                                         |
| <b>Target Names</b>        | GIT1                                                                                                                                                                                                        |

## Image



### Western Blot

Positive WB detected in: Mouse brain tissue


All lanes: GIT1 antibody at 3.5µg/ml

### Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 85, 20, 86 kDa


Observed band size: 85 kDa



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA009439LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA009439LA01HU at dilution of 1:100



Immunofluorescent analysis of Hela cells using CSB-PA009439LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.