



GIT1 Antibody, Biotin conjugated

Product Code	CSB-PA009439LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2X7
Immunogen	Recombinant Human ARF GTPase-activating protein GIT1 protein (470-642AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ARF GTPase-activating protein GIT1 (ARF GAP GIT1) (Cool-associated and tyrosine-phosphorylated protein 1) (CAT-1) (CAT1) (G protein-coupled receptor kinase-interactor 1) (GRK-interacting protein 1), GIT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GIT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.