



GLB1 Antibody, FITC conjugated

Product Code	CSB-PA009476LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16278
Immunogen	Recombinant Human Beta-galactosidase protein (317-455AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-galactosidase (EC 3.2.1.23) (Acid beta-galactosidase) (Lactase) (Elastin receptor 1), GLB1, ELNR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GLB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.