



GNB1 Antibody, FITC conjugated

Product Code	CSB-PA009602LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62873
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-1 protein (157-248AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-1 (Transducin beta chain 1), GNB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GNB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.