



GPR4 Antibody, HRP conjugated

Product Code	CSB-PA009813LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46093
Immunogen	Recombinant Human G-protein coupled receptor 4 protein (285-362AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled receptor 4 (G-protein coupled receptor 19), GPR4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR4