



GPX2 Antibody, HRP conjugated

Product Code	CSB-PA009867LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18283
Immunogen	Recombinant Human Glutathione peroxidase 2 protein (71-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione peroxidase 2 (GPx-2) (GSHPx-2) (EC 1.11.1.9) (Gastrointestinal glutathione peroxidase) (Glutathione peroxidase-gastrointestinal) (GPx-GI) (GSHPx-GI) (Glutathione peroxidase-related protein 2) (GPRP-2), GPX2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GPX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.