



GPX5 Antibody, Biotin conjugated

Product Code	CSB-PA009870LD01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O18994
Immunogen	Recombinant Pig Epididymal secretory glutathione peroxidase protein (22-219AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epididymal secretory glutathione peroxidase (EC 1.11.1.9) (Epididymis-specific glutathione peroxidase-like protein) (EGLP) (Glutathione peroxidase 5) (GPx-5) (GSHPx-5), GPX5
Immunogen Species	Sus scrofa (Pig)
Research Area	Others
Target Names	GPX5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.