



GRHPR Antibody, HRP conjugated

Product Code	CSB-PA009897HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBQ7
Immunogen	Recombinant Human Glyoxylate reductase/hydroxypyruvate reductase protein (1-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyoxylate reductase/hydroxypyruvate reductase (EC 1.1.1.79) (EC 1.1.1.81), GRHPR, GLXR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRHPR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.