



GRHPR Antibody, Biotin conjugated

Product Code	CSB-PA009897HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBQ7
Immunogen	Recombinant Human Glyoxylate reductase/hydroxypyruvate reductase protein (1-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyoxylate reductase/hydroxypyruvate reductase (EC 1.1.1.79) (EC 1.1.1.81), GRHPR, GLXR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRHPR