



GZMK Antibody, Biotin conjugated

Product Code	CSB-PA010084LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49863
Immunogen	Recombinant Human Granzyme K protein (27-264AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granzyme K (EC 3.4.21.-) (Fragmentin-3) (Granzyme-3) (NK-tryptase-2) (NK-Tryp-2), GZMK, TRYP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GZMK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.