



GZMM Antibody, FITC conjugated

Product Code	CSB-PA010085LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51124
Immunogen	Recombinant Human Granzyme M protein (98-243AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granzyme M (EC 3.4.21.-) (Met-1 serine protease) (Hu-Met-1) (Met-ase) (Natural killer cell granular protease), GZMM, MET1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GZMM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.