



HSPB6 Antibody, FITC conjugated

Product Code	CSB-PA010837LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14558
Immunogen	Recombinant Human Heat shock protein beta-6 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock protein beta-6 (HspB6) (Heat shock 20 kDa-like protein p20), HSPB6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPB6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.