




ID4 Antibody

Product Code	CSB-PA010970LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47928
Immunogen	Recombinant Human DNA-binding protein inhibitor ID-4 protein (1-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein inhibitor ID-4 (Class B basic helix-loop-helix protein 27) (bHLHb27) (Inhibitor of DNA binding 4) (Inhibitor of differentiation 4), ID4, BHLHB27
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ID4

Image



Immunofluorescent analysis of HeLa cells using CSB-PA010970LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.