



ID4 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA010970LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P47928 |
| Immunogen | Recombinant Human DNA-binding protein inhibitor ID-4 protein (1-161AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | DNA-binding protein inhibitor ID-4 (Class B basic helix-loop-helix protein 27) (bHLHb27) (Inhibitor of DNA binding 4) (Inhibitor of differentiation 4), ID4, BHLHB27 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | ID4 |