



ID4 Antibody, Biotin conjugated

Product Code	CSB-PA010970LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47928
Immunogen	Recombinant Human DNA-binding protein inhibitor ID-4 protein (1-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein inhibitor ID-4 (Class B basic helix-loop-helix protein 27) (bHLHb27) (Inhibitor of DNA binding 4) (Inhibitor of differentiation 4), ID4, BHLHB27
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ID4