



IDH1 Antibody, Biotin conjugated

Product Code	CSB-PA010989ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75874
Immunogen	Recombinant Human Isocitrate dehydrogenase [NADP] cytoplasmic protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Isocitrate dehydrogenase [NADP] cytoplasmic (IDH) (EC 1.1.1.42) (Cytosolic NADP-isocitrate dehydrogenase) (IDP) (NADP(+)-specific ICDH) (Oxalosuccinate decarboxylase), IDH1, PICD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IDH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.