



IFIT2 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA011021LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P09913 |
| Immunogen | Recombinant Human Interferon-induced protein with tetratricopeptide repeats 2 protein (130-265AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interferon-induced protein with tetratricopeptide repeats 2 (IFIT-2) (ISG-54 K) (Interferon-induced 54 kDa protein) (IFI-54K) (P54), IFIT2, CIG-42 G10P2 IFI54 ISG54 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | IFIT2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |