



IFIT2 Antibody, Biotin conjugated

Product Code	CSB-PA011021LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09913
Immunogen	Recombinant Human Interferon-induced protein with tetratricopeptide repeats 2 protein (130-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced protein with tetratricopeptide repeats 2 (IFIT-2) (ISG-54 K) (Interferon-induced 54 kDa protein) (IFI-54K) (P54), IFIT2, CIG-42 G10P2 IFI54 ISG54
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFIT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.