



IGHD Antibody, Biotin conjugated

Product Code	CSB-PA011113LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01880
Immunogen	Recombinant Human Ig delta chain C region protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunoglobulin heavy constant delta (Ig delta chain C region) (Ig delta chain C region NIG-65) (Ig delta chain C region WAH), IGHD
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IGHD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.