



PRKAG1 Antibody, HRP conjugated

Product Code	CSB-PA01149B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54619
Immunogen	Recombinant Human 5'-AMP-activated protein kinase subunit gamma-1 protein (1-331AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5'-AMP-activated protein kinase subunit gamma-1 (AMPK gamma1) (AMPK subunit gamma-1) (AMPKg), PRKAG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PRKAG1