



IL7R Antibody, Biotin conjugated

Product Code	CSB-PA011670ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16871
Immunogen	Recombinant Human Interleukin-7 receptor subunit alpha protein (21-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-7 receptor subunit alpha (IL-7 receptor subunit alpha) (IL-7R subunit alpha) (IL-7R-alpha) (IL-7RA) (CDw127) (CD antigen CD127), IL7R
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL7R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.