



IP6K1 Antibody, Biotin conjugated

Product Code	CSB-PA011770LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92551
Immunogen	Recombinant Human Inositol hexakisphosphate kinase 1 protein (1-441AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol hexakisphosphate kinase 1 (InsP6 kinase 1) (EC 2.7.4.21) (Inositol hexaphosphate kinase 1), IP6K1, IHPK1 KIAA0263
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IP6K1