



IREB2 Antibody, Biotin conjugated

Product Code	CSB-PA011813LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48200
Immunogen	Recombinant Human Iron-responsive element-binding protein 2 protein (107-216AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iron-responsive element-binding protein 2 (IRE-BP 2) (Iron regulatory protein 2) (IRP2), IREB2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	IREB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.