



KIDINS220 Antibody, HRP conjugated

Product Code	CSB-PA012306LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULH0
Immunogen	Recombinant Human Kinase D-interacting substrate of 220 kDa protein (1077-1370AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinase D-interacting substrate of 220 kDa (Ankyrin repeat-rich membrane-spanning protein), KIDINS220, ARMS KIAA1250
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	KIDINS220
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.