



KIR3DL1 Antibody, HRP conjugated

Product Code	CSB-PA012364LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43629
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 3DL1 protein (81-336AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 3DL1 (CD158 antigen-like family member E) (HLA-BW4-specific inhibitory NK cell receptor) (MHC class I NK cell receptor) (Natural killer-associated transcript 3) (NKAT-3) (p70 natural killer cell receptor clones CL-2/CL-11) (p70 NK receptor CL-2/CL-11) (CD antigen CD158e), KIR3DL1, CD158E NKAT3 NKB1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR3DL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.