



MAP3K13 Antibody, Biotin conjugated

Product Code	CSB-PA013421LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43283
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 13 protein (506-835AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 13 (EC 2.7.11.25) (Leucine zipper-bearing kinase) (Mixed lineage kinase) (MLK), MAP3K13, LZK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K13